



**FERRO®**

NEW

# TIGA MEDICO

Responsible choice

**BTGM2AVL** Standing washbasin mixer



In toilets with public access, it is worth taking advantage of solutions which, in addition to saving water consumption, guarantee the safety of user and allow you to choose new models with medical handles from the Tiga line.



## Safety and sterility



### TIGA MEDICO



#### Standing washbasin mixer

FerroAirMix aerator  
G5/4 drain plug klik-clak  
Ferro Easy Clean System  
chrome

**BTGM2VL**



#### Standing washbasin mixer


swivel spout  
FerroAirMix aerator  
Ferro Easy Clean System  
chrome

**BTGM2AVL**

## Specialist handle


### Safe adjustment

The temperature and flow can be adjusted comfortably and safely with the elbow or wrist and the hands do not come into contact with the mixer surface.



### Ergonomics of the handle

Thanks to the length and design of the handle, the use of the mixer is also easier for people with limited hand efficiency.



There are many advantages of the so-called medical or elbow handle - from those related to hygiene safety to those related to its construction and functionality in use.